

## Anti-Infection Compound Library (96-well)

### Product Details:

Catalog Number:	CS-L002
Formulation:	A collection of 2668 bioactive anti-infection compounds supplied as pre-dissolved Solutions or Solid
Container:	96- or 384-well Plate with Peelable Foil Seal; 96-well Format Sample Storage Tube With Screw Cap and Optional 2D Barcode
Storage:	In solvent: -80°C 2 years, -20°C 1 year; Powder: -20°C 3 years, 4°C 2 years
Shipping:	Blue ice
Packaging:	Inert gas

### Plate layout: CS-L002-1

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Vonoprazan	Peretinoin	AT-130	Bay 41-4109 (racemate)	NSC5844	YYA-021	Faropenem daloxate	Flavopiridol	Flavopiridol (Hydrochloride)	9-Propenyladenine	Empty
b	Empty	(Z)-9-Propenyladenine	Dolutegravir intermediate-1	Tubercidin	Delpazolid	DSM265	E 64c	Aloxistatin	KW-8232	AQ-13 dihydrochloride	Isavuconazonium (sulfate)	Empty
c	Empty	CPI-1189	Nigericin (sodium salt)	Brevianamide F	Cadazolid	Paquinimod	Plerixafor	Aprepitant	AN2718	Nanchangmycin	Dianemycin	Empty
d	Empty	DprE1-IN-2	Golgicide A	Golgicide A-2	Golgicide A-1	NH125	Ferrostatin-1	Isepamicin (sulfate)	Spiramycin	Sodium stibogluconate	Saponins	Empty
e	Empty	Fosfluconazole	Pirlindole	BRD-6929	BM212	4E2RCat	HeE1-2Tyr	Toxoflavin	URB602	Antibiotic-5d	Obefazimod	Empty
f	Empty	CCCP	Cinanserin (hydrochloride)	Flurofamide	Telacebec	Quipazine (dimaleate)	LY294002	Bicyclomycin benzoate	Filibuvir	INCB38579	Sildenafil	Empty
g	Empty	VU0359595	Allopurinol riboside	3-(Methylthio)propionic acid	4-Acetamidobutanoic acid	ST-193	ST-193 (hydrochloride)	PAβN (dihydrochloride)	IT1t (dihydrochloride)	Emodepside	2,4,6-Tribromophenyl caproate	Empty
h	Empty	MBP146-78	PSI-6130	Panidazole	Larixyl acetate	KKL-10	KKL-35	AU1235	KIN1148	Wortmannin	RN-18	Empty

### Plate layout: CS-L002-2

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Antifungal agent 1	PC786	ERDRP-0519	Tanespimycin	Everolimus	Rapamycin	Vorinostat	Ixabepilone	Panobinostat	Telaprevir	Empty
b	Empty	Boceprevir	Danoprevir	Mericitabine	Simeprevir	MK-0608	Narlaprevir	JH-X-119-01	T.cruzi-IN-1	DDD85646	β-Glucuronidase-IN-1	Empty
c	Empty	Enmetazobactam	EHNA (hydrochloride)	Toyocamycin	Reutericyclin	SID 26681509	Trovafoxacin (mesylate)	Fasudil (Hydrochloride)	Takinib	Raltegravir	Raltegravir (potassium)	Empty
d	Empty	GS-441524	M62812	Fumagillol	GW779439X	Canertinib	Canertinib (dihydrochloride)	Gardiquimod	Gardiquimod (diTFA)	Trimetrexate	Sutezolid	Empty
e	Empty	Eperezolid	Linezolid	Emricasan	Ac-CoA Synthase Inhibitor 1	Pocapavir	Remdesivir	AZ960	Balapiravir	R-1479	Terameprocol	Empty
f	Empty	Daclatasvir (dihydrochloride)	Daclatasvir	NM107	Cobicistat	SC75741	Tipifarnib	Iroxacin	T-91825	Misonidazole	Pexiganan	Empty
g	Empty	Rifalazil	Fantofarone	LEDGIN6	AzddMeC	Sulopenem	Betulinic acid	PBT 1033	Tegobuvir	Taribavirin (hydrochloride)	Nevirapine	Empty
h	Empty	Delavirdine (mesylate)	Efavirenz	Rilpivirine	Rilpivirine (hydrochloride)	PD 169316	Gatifloxacin	Gatifloxacin (hydrochloride)	5-Azacytidine	Lexithromycin	Dicresulene (diammonium)	Empty

### Plate layout: CS-L002-3

1	2	3	4	5	6	7	8	9	10	11	12
---	---	---	---	---	---	---	---	---	----	----	----

a	Empty	MMV390048	Vapendavir (diphosphate)	WIN 54954	Duramycin	Zabofloxacin (hydrochloride)	Neticonazole	Eberconazole nitrate	Ceftobiprole medocartil (sodium)	Prussian blue insoluble	Mitoguazone	Empty
b	Empty	Chloroquinoline sulfonamide	Polymyxin B nonapeptide	Polymyxin B nonapeptide (TFA)	3'-Azido-2',3'-dideoxyuridine	Letrazuril	Peldesine (dihydrochloride)	Tetroxoprim	TMC310911	Bacitracin	N-Butylthiophosphoric	Empty
c	Empty	$\beta$ -Cyclodextrin	Polyinosinic-polycytidylic acid	Poly (I:C):Kanamycin (1:1)	Selamectin	Cefathiamidine	(6R,7S)-Cefminox (sodium)	$\beta$ -Chloro-L-alanine	U18666A	OSS_128167	HIV-1 integrase inhibitor 8	Empty
d	Empty	Nosiheptide	Imidocarb (dipropionate)	Conessine	R 59-022 (hydrochloride)	TCS PrP Inhibitor 13	Duocarmycin TM	Fluensulfone	Ganfeborole (hydrochloride)	Ectoine	Potassium guaiaacolsulfonate	Empty
e	Empty	Inosine pranobex	Nicarbazine	A40926	Decamethoxine	Rezafungin (acetate)	ME1111	Ganaplacide (hydrochloride)	Valopicitabine (dihydrochloride)	Z-LVG-CHN2	Bisindolylmaleimide IV	Empty
f	Empty	Tomoglovir	Subactam pivoxil	Miconomicin (sulfate)	T-00127_HEV1	Brequinar	6-Diazo-5-oxo-L-nor-Leucine	Cruzain-IN-1	Gamithromycin	Cefodizime	Cefodizime (sodium)	Empty
g	Empty	exo-IWR-1	Pretomanid	Polygodial	Delamanid	ML-SA1	Loxoribine	TMPyP4 tosylate	PB28 (dihydrochloride)	Arterolane	Alexidine (dihydrochloride)	Empty
h	Empty	HLI373 (dihydrochloride)	PPNDS (tetrasodium)	RO-9187	XZ426	Clotrimazole	Carindacillin (sodium)	Troleandomycin	(Rac)-Modipafant	Nitrocefin	Morphothiadin	Empty

Plate layout: CS-L002-4

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	SDZ285428	GlyRS-IN-1	Demecycline	JFD00244	Afabcin	Enzaplatovir	(S)-Enzaplatovir	Nacubactam	Tenofovir exalidex	Baloxavir marboxil	Empty
b	Empty	Baloxavir	Quilseconazole	Elsulfavirine	Selgantolimo d	Ziresovir	Bersacapavir	Vebicovir	Vonafexor	HZ-1157	BM635	Empty
c	Empty	NITD-349	Vadimezan	Obatoclox (Mesylate)	AN0128	Tavorole	ISPA-28	Terreic acid	(Z)-Guggulsterone	GNF-2	Furamide (dihydrochloride)	Empty
d	Empty	(E/Z)-Teriflunomide	UCM05	TMC353121	Icerguastat	Nilofabacin	MS417	AR-9281	Pyocyanin	TLR7 agonist 1	Tafenoquine (Succinate)	Empty
e	Empty	(3R,4R)-A2-32-01	(3S,4S)-A2-32-01	3'-Azido-3'-deoxy-5-methylcytidine	3'-Azido-3'-deoxy-5-fluorocytidine	6-Amino-5-azacytidine	CWHM-1008	TBAJ-587	Jaspamycin	3,4'-Dihydroxyflavone	ACT-451840	Empty
f	Empty	MCB-3681	ABBV-744	GSK2200150 A	AB-423	Rilematovir	ZINC0312931 g	Terizidone	UCT943	S119-8	JNJ-632	Empty
g	Empty	MBX-4132	BRITE-338733	Dextran	Aristeromycin	RSV-IN-1	LHF-535	(E)-LHF-535	Fosmidomycin sodium salt	Gln-AMS (TFA)	Arg-AMS	Empty
h	Empty	Prostaglandin A2	Glycolithocholic acid 3-sulfate	Hydroxyphenyl lactic acid	Brassicasterol	5'-Deoxyadenosine	Glycogen, Mussel	BC-7013	ASN0357680 0	N-Butanoyl-DL-homoserine	DSHS00884	Empty

Plate layout: CS-L002-5

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Actinonin	Tuberculosis inhibitor 3	T-2307	BA-53038B	AWZ1066S	N-3-oxo-dodecanoyl-L-homoserine	DB772 (hydrate)	Retro-2 cycl	4-(tert-Butyl)-benzhydroxamic Acid	Azamethipho s	Empty
b	Empty	Metronidazole acetic acid	Parbendazole	CRS3123 (dihydrochloride)	Bio-AMS	BRD5529	Lufenuron	VU0420373	Bisdionin C	VP-4604	Oleandomycin	Empty
c	Empty	Conoidin A	Picaridin	Omaciclovir	Sparstolonin B	Cyprodinil	Cadroxifloxacin	Dextran sulfate sodium salt	Dextran sulfate sodium salt	Dextran sulfate sodium salt	Metaflumizone	Empty
d	Empty	Viridicatinol	Lotilaner	Halocarban	116-9e	BLT-1	Lipofermata	Lalistat 1	MAC13772	5-Pentadecylresorcinol	TH1020	Empty
e	Empty	N-(Hydroxymethyl)nicotinamide	IR415	Purfalcamine	GRL0617	Tetraconazole	Diallyl Trisulfide	NGI-1	Methoxyfenozide	Coblopassvir (dihydrochloride)	LP-935509	Empty
f	Empty	RG7834	Lincomycin	L-(+)-threo-Chloramphenicol	Cabamiquine	Cabamiquine (succinate)	IHVR-11029	PC945	Fenhexamid	(R)-ZINC-3573	Florfenicol amine (hydrochloride)	Empty
g	Empty	Fialuridine	PKR-IN-C51	EcDsbB-IN-9	2,4-Diacetylphloroglucinol	VU041	Benzalpthalide	Anhydrotetracycline (hydrochloride)	N-dodecanoyl-L-Homoserine	MreB Perturbing Compound	SNX-0723	Empty
h	Empty	GSK983	K777	Fluopyram	Camalexin	Propineb	Epoxiconazole	Bifenazate	Fosmanogepix	Trofosamide	2-Hydroxychalcone	Empty

Plate layout: CS-L002-6

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Diloxanide	Boscalid	Cyclopyrimor ate	GNF6702	PF-3450074	R-10015	Iniparib	U0126-ETOH	RYL-552	PIK-93	Empty
b	Empty	ASP6432	N- Acetyltyramin e	Hesperadin	M4284	Sabizabulin	BMS-433771	BMS-433771 (dihydrochlori de hydrate)	Quinoxyfen	GSK369796 Dihydrochlori de	Zidebactam	Empty
c	Empty	Chalcone	Gloxazone	Alizarin complexone	1,8- Dihydroxy- 4,5-	Cefazedone	Dehydroemeti ne	Demeclocycli ne	Roseoflavin	L-Ribulose	Brodimoprim	Empty
d	Empty	Evernic Acid	Resorantel	BKI-1369	3- Methoxybenz amide	Torcitabine	Amylmetacre sol	α-Thujone	Monolaurin	Inz-5	DDX3-IN-1	Empty
e	Empty	DDX3-IN-2	Cyclo(L-Phe- L-Val)	ML-60218	IMB-301	Clovamide	Kyotorphin	Umrirolimus	Stampidine	Retro-2	Aurintricarbox ylic acid	Empty
f	Empty	AVG-233	Surfen (dihydrochlori de)	FR 180204	bpV(phen) (trihydrate)	Rostafuroxin	Metronidazole Benzoate	Triadimefon	Q-VD-OPH	Pyrimorph	Trifloxystrobin	Empty
g	Empty	HDAC3-IN- T247	5- Hydroxymebe ndazole	Sisunatovir	Sisunatovir (hydrochlorid e)	H1PVAT	Pimodivir	MMV008138	Quinocetone	PF 1022A	VT-1598	Empty
h	Empty	Aminothiazole	Diminazene (aceturate)	K6PC-5	Dodecyl gallate	N-Octanoyl-L- homoserine lactone	Slingshot inhibitor D3	TPCK	CA inhibitor 1	NVR 3-778	GLP-26	Empty

Plate layout: CS-L002-7

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	FGI-106 (tetrahydrochl oride)	GPI-1046	G0507	IHVR-19029	ML406	TBT1	GAK inhibitor 49	GAK inhibitor 49 (hydrochlorid e)	Qstatin	Epetraborole (hydrochlorid e)	Empty
b	Empty	ABMA	TTP-8307	TBA-354	LASV inhibitor 3.3	iKIX1	Hck-IN-1	GSK-A1	G907	Oltipraz	LEI110	Empty
c	Empty	Velpatasvir	Rifamycin S	2'-C- Methyladenos ine	WCK-4234	Olanexidine (hydrochlorid e)	Ferutinin	Kresoxim- methyl	2- Mercaptopyri dine N-oxide	2,3-Dehydro- 2-deoxy-N- acetylneurami n	Cytidine-5'- triphosphate	Empty
d	Empty	Vulpinic acid	Paritaprevir	PNU-101603	ML251	(±)-Alliin	KIN101	Desacetylcefo taxime	LysRs-IN-2	Acetylvostat in	NPD-1335	Empty
e	Empty	TCMDC- 135051	TCMDC- 135051 (hydrochlorid e)	Dichlorophen	Bephenium (hydroxynaph thoate)	Pyrantel (pamoate)	Pyrantel (tartrate)	Diethylcarba mazine (citrate)	ZL0580	Eprinomectin	Myristoyl glutamic acid (sodium)	Empty
f	Empty	FASN-IN-4 (tosylate)	Phleomycin	Primaquine (diphosphate)	Primaquine	LDC4297	LDC4297 (hydrochlorid e)	Piperlongumi nine	NP213 (TFA)	Tizoxanide	HBF-0259	Empty
g	Empty	GSK3494245	Guanabenz (hydrochlorid e)	ML324	Trimethobenz amide (hydrochlorid e)	6-OAU	Cycloguanil (hydrochlorid e)	Albendazole sulfoxide	5-(N,N- Hexamethyl ene)-amiloride	Sibofimloc	AN3661	Empty
h	Empty	AEBSF (hydrochlorid e)	RNPA1000	IMB-XH1	KH-CB19	AN7973	Ibezapolstat	Neticonazole (hydrochlorid e)	Benzylidodecy ldimethylamm onium	2,3- Butanediol	Tylvalosin (tartrate)	Empty

Plate layout: CS-L002-8

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Allyl methyl sulfide	Dehydroaceti c acid (sodium)	Fusaric acid	N-Desethyl amodiaquine	N-Desethyl amodiaquine dihydrochlorid e	Phosphonoac etic acid	MRL-494 (hydrochlorid e)	hDHODH-IN- 4	N-Acetyl-D- mannosamin e	DNA31	Empty
b	Empty	Cefminox (sodium)	N- Glycolylneuram inic acid	LHVS	Macozinone	Ramoplanin	TCA1	3-Formyl rifamycin	Amoscanate	Isoeulertherin	Flutrimazole	Empty
c	Empty	Nourseothrici n (sulfate)	FR179642	EACC	KDU691	ML318	Stearyl-dietha nolamine	SPR719	N- (Ketocaproyl)- DL-	4α,7α,7αc- Nepetalacton e	PTACH	Empty
d	Empty	Surfactin	Asperphenam ate	Ferulenol	Pimozide	RSV604	RSV604 (racemate)	Benzoethiohyd razide	5- Hydroxymeth yl-2'-	Maraviroc	HQNO	Empty
e	Empty	Sapientic acid (sodium)	Reverse transcriptase- IN-1	HIV-1 integrase inhibitor	Depufavirine	MAC-545496	Daunorubicin (hydrochlorid e)	LpxH-IN-AZ1	JH-LPH-28	JH-LPH-33	LoICDE-IN-1	Empty
f	Empty	Furanone C- 30	HADA (hydrochlorid e)	CBS1117	MBX2329	Mequindox	AMOZ	Ristomycin (sulfate)	Anti-virus agent 1	Hydroxychlor oquine Impurity E	HBV-IN-4	Empty
g	Empty	LXE408	Emtricitabine triphosphate (tetrasodium)	FNDR-20123	Dynapyrazole -A	Diethylamine NONOate (diethylammo nium)	BVDV-IN-1	TP0586532	Desacetylcep haporin (sodium)	ONX-0914	ONX-0914 (TFA)	Empty
h	Empty	Zanamivir	HCoV-229E- IN-1	DSM502	DSM705 (hydrochlorid e)	CHIKV-IN-2	YM-201636	Thiamphenico l glycinate (hydrochlorid e)	HIV-1 inhibitor-8	CCF0058981	CDK9-IN-1	Empty

Plate layout: CS-L002-9

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Rifaximin	Dolutegravir	Dolutegravir (sodium)	BMS-707035	SARS-CoV-2-IN-6	Avibactam tomilopil	Pseudomonic acid C	Azadirachtin B	PK150	Osetamivir	Empty
b	Empty	Osetamivir acid	Anguizole	BMS-986094	N-Hexanoyl-L-homoserine lactone	Chamaechromone	GRP-60367 (hydrochloride)	Levofloxacin N-oxide	TAK-779	Isoflupredone acetate	Tubacin	Empty
c	Empty	Ionomycin (calcium)	PF429242 (dihydrochloride)	Finafloxacin	SM-102	VCH-916	D75-4590	NITD-2	GST-HG131	Cyclotriazadisulfonamide	HIV-1 inhibitor-6	Empty
d	Empty	SID-852843	Quabodepistat	Mitoxantrone	Mitoxantrone (dihydrochloride)	N,N',N"-Triacetylchitrose	Glycerol monocaprate	Camostat (mesylate)	Aloperine	(-)-Cedrene	Cefcapene pivoxil (hydrochloride)	Empty
e	Empty	AG-1478	(+)-Neomenthol	Amphotericin B methyl ester	Desmethyl Levofloxacin	3-Desacetyl Cefotaxime lactone	(R)-?Linezolid	Decarboxy Moxifloxacin	Nifurpirinol	Anidulafungin	Fluxapyroxad	Empty
f	Empty	Ascomycin	hDHODH-IN-3	Artelinic acid	AVN-944	PFDHODH-IN-1	Hexyl gallate	BSH-IN-1	Biapenem	Penconazole	Brivudine	Empty
g	Empty	BTZ043	BTZ043 (Racemate)	Cletoquine (oxalate)	Desethyl chloroquine (diphosphate)	LtaS-IN-1	Carbendazim	Aspoxicillin	Desmethyl ferroquine	Molnupiravir	SARS-CoV-IN-1	Empty
h	Empty	SARS-CoV-IN-2	SARS-CoV-IN-3	Cefsulodin (sodium)	TH1217	Chlorin e6	Cytarabine	Cytarabine (hydrochloride)	Mefentrifluconazole	Mpro inhibitor N3 (hemihydrate)	Methyl carnosate	Empty

Plate layout: CS-L002-10

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	SIRT1-IN-1	IA-Alkyne	Entecavir	Entecavir (monohydrate)	Epirubicin (hydrochloride)	Ertapenem sodium	BRD0539	HCV-IN-30	XP-59	Etoposide	Empty
b	Empty	X77	Etoposide phosphate	Metallo $\beta$ -lactamase ligand 1	Oxazosulfil	Picoxystrobin	Ganciclovir	Ganciclovir (sodium)	N-Decanoyl-L-homoserine lactone	Ethylhydrocupreine (hydrochloride)	Toltrazuril sulfoxide	Empty
c	Empty	Hydroxymetronidazole	4-Epitetracycline	CYM50358	A2ti-1	A2ti-2	Coenzyme FO	6-Azathymine	2-Hydroxy-4-methylbenzenesulphonic	Remdesivir O-desphosphat	$\omega$ -Hydroxy-DEET	Empty
d	Empty	Lansoprazole	NPD-001	Levamisole (hydrochloride)	LY2334737	Megestrol acetate	Meropenem	Meropenem (trihydrate)	Miltefosine	Nilutamide	PF-00835231	Empty
e	Empty	Oglufanide	Pirarubicin (Hydrochloride)	Neuraminidase-IN-1	6'-Sialyllactose (sodium)	Quinacrine (dihydrochloride)	l-Atabrine dihydrochloride	d-Atabrine dihydrochloride	Lactonic sophorolipid	Resiquimod	Fanotaprim	Empty
f	Empty	Modoflaner	EBOV/MARV-IN-1	Ebselen	Tacrolimus (monohydrate)	DMT-dl Phosphoramidite	Tetrandrine	6-Thioguanine	Valrubicin	Zoledronic Acid	Tenofovir Disoproxil (fumarate)	Empty
g	Empty	Tenofovir Disoproxil	Pirodavir	GRL-0496	Bemnifosbuvir (hemisulfate)	Bemnifosbuvir	Questin	Fexinidazole	SSAA09E2	Lufotrelvir	SSAA09E3	Empty
h	Empty	Ticlatone	ALC-0315	$\beta$ -Lactamase-IN-2	HCV-IN-31	Atovaquone	ELQ-300	Evixapodlin	1-Dodecylimidazole	2'-O-Me-C(Bz) Phosphoramidite	2'-OMe-G(ibu) Phosphoramidite	Empty

Plate layout: CS-L002-11

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	2'-OMe-A(Bz) Phosphoramidite	DMT-dA(PAc) Phosphoramidite	DMT-dG(dmF) Phosphoramidite	DMT-dC(ac) Phosphoramidite	7-TFA-ap-7-Deaza-dG	Clevudine	5'-O-TBDMS-N4-Benzoyl-2'	5'-O-TBDMS-Bz-dA	5'-O-TBDMS-dT	5'-O-TBDMS-dG	Empty
b	Empty	5'-O-TBDMS-dA	5'-O-DMT-N4-Bz-5-Me-dC	5'-O-DMT-N4-Ac-2'-F-dC	5'-O-DMT-N4-Bz-2'-F-dC	5'-O-DMT-N6-Me-2'-dA	5'-O-DMT-N2-DMF-dG	5'-O-DMT-rU	5'-O-DMT-Bz-rC	5'-O-DMT-2'-O-TBDMS-Bz-rC	5'-O-DMT-3'-O-TBDMS-Ac-rC	Empty
c	Empty	5'-O-DMT-2'-O-TBDMS-rI	5'-O-DMT-2'-O-TBDMS-Ac-rC	Antibacterial agent 19	Monocaprylin	Nirmatrelvir	5-Benzoyloxygramine	Octaethylene glycol monododecyl	Merafloxacin	SU0268	IL-4-inhibitor-1	Empty
d	Empty	Tenofovir	Tenofovir (hydrate)	iHCK-37	YH-53	TBI-223	SQ609	RdRP-IN-2	TL-895	SRPIN340	Zapnometinib	Empty
e	Empty	Temocillin (disodium)	(+)-JNJ-A07	Dup-721	ZCL278	Antiviral agent 5	Merimepodib	MERS-CoV-IN-1	Ombitasvir	Dasabuvir	SDH-IN-1	Empty
f	Empty	BMS-378806	Rimonabant	Rimonabant (Hydrochloride)	D-(-)-3-Phosphoglyceric acid	PI-55	Cotrimoxazole (Trimethoprim)	SDZ 224-015	BRD-8000.3	Antibacterial agent 26	Antibacterial agent 27	Empty
g	Empty	SARS-CoV-2-IN-7	DC-LC3in-D5	Z-Leu-Arg-AMC	Cymal-6	Dapivirine	Lersivirine	Ravuconazole	Isavuconazole	Lanoconazole	(Z)-Lanoconazole	Empty

h	Empty	Luliconazole	Cimetidine	Ensitrelvir	Ensitrelvir (fumarate)	FBA-IN-1	Emodin	Indomethacin	Indomethacin (sodium hydrate)	CCR8 antagonist 1	Cipargamin	Empty
---	-------	--------------	------------	-------------	------------------------	----------	--------	--------------	-------------------------------	-------------------	------------	-------

**Plate layout: CS-L002-12**

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Asunaprevir	DNMT3A-IN-1	MAC13243	FPR2 agonist 2	NS2B/NS3-IN-3	PA3552-IN-1	SARS-CoV-2-IN-13	SARS-CoV-2-IN-14	Bunamidine (hydrochloride)	AAA-10	Empty
b	Empty	Methotrexate	Methotrexate (disodium)	Flupyradifurone	Triperiden	Salazopyridazine	Aripiprazole	Claficapavir	Linvecorvir	Lopinavir	Kaempferol	Empty
c	Empty	DMA-135 (hydrochloride)	INE963	Ribavirin carboxylic acid	Obeldesivir	Clioquinol	EIDD-2749	PqsR/LasR-IN-2	PqsR/LasR-IN-3	Mtb ATP synthase-IN-1	Dexamethasone	Empty
d	Empty	Dexamethasone acetate	Aspirin	Sulfasalazine	Thalidomide	SARS-CoV-2 3CLpro-IN-2	GSK2556286	Dieckol	PRRSV/CD163-IN-1	LpxC-IN-10	Acloproxalap	Empty
e	Empty	Dovramilast	HIV-IN-6	Osetamivir acid methyl ester	ANT3310 (sodium)	Tribendimidine	Ceftaroline fosamil	Elvitegravir	Ulonivirine	Vipoglanstat	Golotimod	Empty
f	Empty	Golotimod (hydrochloride)	Pyronaridine (tetraphosphate)	Favipiravir	Nesbuvir	Levomefolic acid	Radezolid	Amenamevir	ACSS2-IN-2	Delafloxacin	Delafloxacin (meglumine)	Empty
g	Empty	Diamthazole (hydrochloride)	Laninamivir	NLRP3-IN-11	Laninamivir octanoate	Gentamicin C1a	GLUT1-IN-2	IRG1-IN-1	AB-836	UCK2 Inhibitor-1	Antimicrobial agent-14	Empty
h	Empty	MMV688845	Tedizolid	Tedizolid (phosphate)	ccc_R08	TBI-166	CamA-IN-1	1-Deoxyynojirimycin	Omadacycline (tosylate)	Omadacycline (hydrochloride)	Avibactam (free acid)	Empty

**Plate layout: CS-L002-13**

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Rolusafine	Avibactam (sodium)	Avibactam (sodium hydrate)	Bedaquiline	Bedaquiline (fumarate)	(Rac)-Bedaquiline	Cenicriviroc	Cenicriviroc Mesylate	AT-533	D-Aminophenyl mercuric	Empty
b	Empty	Umifenovir (hydrochloride)	Vidofludimus	Dexelvucitabine	Levonadifloxacin	Pafuramidine	Ozenoxacin	AX20017	SQ109	Sofosbuvir	ClpB-IN-1	Empty
c	Empty	TH-Z93	TH-Z145	LabMol-319	Sildenafil	Sildenafil (citrate)	ARN1468	Hispolon	SARS-CoV-2 nsp13-IN-1	Mab?Aspartate Decarboxylase	Chitin synthase inhibitor 4	Empty
d	Empty	NDM-1 inhibitor-3	HIV-1 inhibitor-45	XR8-89	ATP synthase inhibitor 2	3-Chlorogentisyl alcohol	Antitubercular agent-30	Lonafarnib	Staurosporine	Doxorubicin (hydrochloride)	PIKfyve-IN-1	Empty
e	Empty	Tipranavir	Elaiophylin	Antibacterial agent 125	GSK 650394	Antiviral agent 24	Acetyl-CoA Carboxylase-IN-1	CLK1-IN-2	CB2R/FAAH modulator-3	Geldanamycin	Tenofovir alafenamide	Empty
f	Empty	Tenofovir alafenamide fumarate	Tenofovir alafenamide hemifumarate	Letermovir	PSI-6206	MenA-IN-2	UNC0638	BMS-265246	Epothilone D	E-64	Nelfinavir	Empty
g	Empty	Nelfinavir (Mesylate)	Vesnarinone	Grazoprevir	Grazoprevir potassium salt	Grazoprevir hydrate	MMV676584	Pritelivir	Pritelivir (mesylate)	Pritelivir (mesylate hydrate)	Dynasore	Empty
h	Empty	Eltrombopag	Eltrombopag (Olamine)	Flaviviruses-IN-3	Avermectin B1a	Transketolase-IN-4	Ivermectin	Avermectin B1	Enpp-1-IN-16	AG 555	WU-FA-01	Empty

**Plate layout: CS-L002-14**

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Trovirdine	Emivirine	Loviride	BIO-acetoxime	ARC7	MMV009085	Sacubitril	Temsavir	Fostemsavir	Fostemsavir Tris	Empty
b	Empty	RDR 02308	Kaempferide	Roflumilast	Tricribine	CHIR-090	Imatinib	Salubrinal	Scriptaid	Tyrphostin A9	Braco-19	Empty
c	Empty	Braco-19 (trihydrochloride)	Tilbroquinol	TMC647055 (Choline salt)	Cabotegravir	Cabotegravir (sodium)	Salinomycin	Vesatolimod	Ledipasvir	Ledipasvir (acetone)	Ledipasvir (D-tartrate)	Empty
d	Empty	Efinaconazole	Tulathromycin A	Puromycin aminonucleoside	Ansamitocin P-3	Morinidazole	(R)-Morinidazole	Elbasvir	Bromosporine	FTI-277 (hydrochloride)	GSK2838232	Empty
e	Empty	AMD 3465 (hexahydrobromide)	GNF179	Cyclacillin	Maribavir	Micafungin (sodium)	Nimorazole	RRx-001	Sulfacytine	Camptothecin	Resveratrol	Empty
f	Empty	cis-Resveratrol	Deferasirox (Fe3+ chelate)	Kanamycin (sulfate)	Oligomycin A	Birinapant	Brefeldin A	Apabetalone	SN 2	ESI-09	Dapaconazole	Empty

g	Empty	Sarolaner	Gepotidacin	Relebactam	Ridnirilazole	Artefenomel	Avarofloxacin	Doravirine	Censavudine	(2S,5S)-Censavudine	Fosravuconazole	Empty
h	Empty	Fosravuconazole (L-lysine ethanolate)	Samatasvir	Rolipram	Lefamulin (acetate)	AFN-1252	CK-666	5-Methylthioadenosine	LJ001	PD 404182	Fluralaner	Empty

Plate layout: CS-L002-15

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Afoxolaner	Saquinavir (mesylate)	Saquinavir	Toltrazuril (sulfone)	Retapamulin	Peramivir (trihydrate)	Oseltamivir (phosphate)	Rifabutin	Epothilone B	Doramectin	Empty
b	Empty	Naphthoquine (phosphate)	Darunavir	Darunavir (Ethanolate)	Loratadine	Nepafenac	Loteprednol Etabonate	Deferasirox	Vancomycin (hydrochloride)	Dimethyl fumarate	Atazanavir (sulfate)	Empty
c	Empty	Posaconazole	Vicriviroc (maleate)	Idarubicin (hydrochloride)	Loxapine	Loxapine (succinate)	Zalcitabine	Terbinafine hydrochloride	Terbinafine	Butenafine (Hydrochloride)	Mevastatin	Empty
d	Empty	Limonin	Minocycline (hydrochloride)	Zidovudine	Tenatoprazole	Acyclovir	Abacavir	Abacavir (sulfate)	Penciclovir	Valacyclovir (hydrochloride)	Famciclovir	Empty
e	Empty	Emtricitabine	Amprenavir	Fosamprenavir (Calcium Salt)	Mefloquine (hydrochloride)	Cidofovir	Salinomycin (sodium salt)	Azathramycin	Sivelestat	Sivelestat (sodium)	Sivelestat (sodium tetrahydrate)	Empty
f	Empty	Aplaviroc (hydrochloride)	Cefditoren (Pivoxil)	Garenoxacin (Mesylate hydrate)	Sancycline	Sancycline (hydrochloride)	Mizoribine	Metoprolol	Rosuvastatin (Calcium)	Rosuvastatin (Sodium)	Azithromycin	Empty
g	Empty	Pantoprazole (sodium)	Pantoprazole (sodium hydrate)	Clarithromycin	Itraconazole	Tolfenpyrad	Valifenalate	(S,S)-Valifenalate	Penthiopyrad	Meptyldinocap	Haloxon	Empty
h	Empty	PLpro inhibitor	Demeclocycline (hydrochloride)	G-418 (disulfate)	Capreomycin (sulfate)	Heparin (Lithium salt)	Thio-TEPA	Berberine (chloride hydrate)	Pneumocandin B0	Fidaxomicin	Buparvaquone	Empty

Plate layout: CS-L002-16

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Griseofulvin	Dalbavancin (hydrochloride)	Chloroquine (phosphate)	Chloroquine	Bithionol	Bithionol (sulfoxide)	Solithromycin	Oxyclozanide	Mebendazole	Closetel	Empty
b	Empty	Closetel (sodium)	Febantel	Rafoxanide	Piperazine (citrate)	Bictegravir	Bictegravir (sodium)	Framycetin (sulfate)	Glecaprevir	Oteseconazole	Zoliflodacin	Empty
c	Empty	Ochromycino	Pyrimethamine	MUT056399	Walrycin B	Leupeptin (hemisulfate)	Rolitetracycline	Berberine (chloride)	Purvalanol B	IOWH-032	4-Phenoxybenzylamine	Empty
d	Empty	Chitinase-IN-1	(±)-BI-D	PFI-2 (hydrochloride)	Galidesivir (hydrochloride)	Galidesivir	SIBA	Targocil	BRD7116	Ornidazole (Levo-)	Dextrorotation nimorazole phosphate	Empty
e	Empty	Endoxifen (hydrochloride)	BQR-695	FIT-039	Rotterlin	Anisomycin	Ro24-7429	Sulfaclozine	Sulfaclozine (sodium)	Azvudine	Azvudine (hydrochloride)	Empty
f	Empty	OG-L002	Ferroquine	Besifovir	Ribocil	Ribocil-C	Artemisone	Ladarixin (sodium)	PK 11195	Teclozan	BRD3308	Empty
g	Empty	STING agonist-1	Miransertib	FOY 251	Triazavirin	LED209	TBA-7371	β-Lactamase-IN-1	Voxilaprevir	Brilacidin (tetrahydrochloride)	ITX5061	Empty
h	Empty	Acoziborole	Contezolid	Vaborbactam	ACHN-975 (TFA)	Pleconaril	KIN1408	Palmitoylethanolamide	2-Phenylethyl isothiocyanate	D-?(+)-?Phenylactic acid	Clemizole (hydrochloride)	Empty

Plate layout: CS-L002-17

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Benzylamine (hydrochloride)	Phenazine-1-carboxylic acid	2'-Hydroxy-4'-methylacetophenone	Tofacitinib (citrate)	Ganciclovir mono-O-acetate	2,3-Dimethoxybenzaldehyde	5'-DMT-3'-TBDMS-ibuprofen	Oseltamivir-acetate	Diethyl butylmalonate	S.pombe lumazine synthase-IN-1	Empty
b	Empty	Genevant CL1	CP26	Nucleozin	Mavoxiafor (trihydrochloride)	Asparagusic acid	Temsirolimus	Plerixafor (octahydrochloride)	4-Acetamidophenyl acetate	D609	SPK-601	Empty
c	Empty	Lomibuvir	Voriconazole	5-Hydroxypyrazine-2-	Faropenem sodium	I2906	NBD-556	NBD-557	D-Ribonolactone	Furagin	Demethoxynececalin	Empty
d	Empty	Nucleoside-Analog-1	Nucleoside-Analog-2	Benzyl isothiocyanate	2-Methoxybenzaldehyde	Ibuprofen	Ibuprofen (sodium)	Fosamprenavir	S-Methylisothiourea (sulfate)	(Tetrahydro-2H-pyran-4-yl)methanol	Methyl indole-3-carboxylate	Empty
e	Empty	Ritonavir	Etravirine	5-Fluorouracil	Pentostatin	Valganciclovir (hydrochloride)	Nifuratel	Trifluridine	Telithromycin	Aurothioglucose	Tildipirosin	Empty



f	Empty	Tebipenem	Netilmicin (sulfate)	Cefotaxime (sodium)	Cefotaxime	Colistin (sulfate)	Nitrofurantoin	Thiamine monochloride	Levamisole	Cefotiam hexetil (hydrochlorid)	Cefetamet	Empty
g	Empty	Plicamycin	Sulfalene	Chlorphenesi n	Halofantrine (hydrochlorid e)	Cephapirin	Cephapirin (sodium)	Quinupristin	Cilastatin	Cilastatin (sodium)	Trovafloracin	Empty
h	Empty	Adenosine 5'- monophosphate	Rosoxacin	Cerulenin	Colistin methanesulfo nate (sodium)	Flucloxacillin sodium	Polymyxin B (Sulfate)	Cefpodoxime	Cefacetile	Clavulanic acid	Clavulanate (lithium)	Empty

Plate layout: CS-L002-18

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Mecillinam	Gentamicin (sulfate)	Cyproconazole	Hexaconazole	Pristinamycin	Pristinamycin IA	Ertapenem (disodium)	Toremifene (citrate)	Carvedilol	Pamidronic acid	Empty
b	Empty	Telbivudine	Prulifloxacin	Ceftiofur (hydrochlorid e)	Valnemulin (hydrochlorid e)	D-Cycloserine	Sulfamethazine	Azulene	Leflunomide	Artemisinin	Floxuridine	Empty
c	Empty	Fluconazole	Ketoconazole	(+)- Ketoconazole	(-)- Ketoconazole	Daptomycin	Omeprazole	Omeprazole (sodium)	Stavudine	Tigecycline	Tigecycline (tetramesylat e)	Empty
d	Empty	Ofloxacin	Marbofloxacin	Aztreonam	Norfloxacin	Natamycin	Bestatin	Bestatin (hydrochlorid e)	Cefdinir	Flucytosine	Estradiol	Empty
e	Empty	Alpha-Estradiol	Pefloxacin (mesylate)	Nicotinamide	Fluticasone (propionate)	Ketotifen (fumarate)	Olsalazine (Disodium)	Toltrazuril	Tinidazole	Imiquimod	Imiquimod (hydrochlorid e)	Empty
f	Empty	Imiquimod (maleate)	Carmofur	Cefoselis (sulfate)	Doripenem (monohydrate )	Nafamostat (mesylate)	Cefaclor	Cefaclor (monohydrate )	Mycophenolate Mofetil	Cephalexin	Cephalexin (monohydrate )	Empty
g	Empty	Cefoperazone	Cefoperazone (sodium salt)	Sulfapyridine	Sulfameter	Acetylcysteine	Nitazoxanide	Allopurinol	Erythromycin	Erythromycin (thiocyanate)	Amphotericin B	Empty
h	Empty	1-Docosanol	Albendazole	Nitrofurazone	Amorolfine (hydrochlorid e)	Sulfanilamide	Praziquantel	Didanosine	Lamivudine	Adefovir dipivoxil	Methylprednisolone	Empty

Plate layout: CS-L002-19

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Thiabendazole	Enoxacin (hydrate)	Rifapentine	Pyrazinamide	Rifampicin	Sulfadiazine	Sulfadiazine (sodium)	Chlorprothixene	Chlorprothixene (hydrochlorid)	Oxytetracycline	Empty
b	Empty	Oxytetracycline (hydrochlorid)	Ethionamide	Vidarabine	Vidarabine phosphate	Ranitidine (hydrochlorid e)	Erdosteine	Oxfendazole	Butoconazole (nitrate)	Flubendazole	Chloroxine	Empty
c	Empty	Oxibendazole	Bifonazole	Diphenhydramine (hydrochlorid)	Prothionamide	Idoxuridine	Sparfloxacin	Hydroxyurea	Metronidazole	Tioconazole	Sulfamethoxazole	Empty
d	Empty	Sulfisoxazole	Triamcinolone	Isoniazid	Levofloxacin	Levofloxacin (hydrate)	Levofloxacin (hydrochlorid e)	Sulfamethizole	Sulbactam	Sulbactam (sodium)	Sulfadimethoxine	Empty
e	Empty	Sulfadimethoxine (sodium)	Rimantadine (hydrochlorid e)	Sarafloxacin (hydrochlorid e)	Liranaftate	Ciprofloxacin	Ciprofloxacin (hydrochlorid e)	Diclazuril	Fenticonazole (Nitrate)	Dyclonine (hydrochlorid e)	Tolnaftate	Empty
f	Empty	Bromhexine (hydrochlorid e)	Mecarbinat	Probuco	Sitafoxacin (hydrate)	Tebipenem pivoxil	Nalidixic acid	D-Sorbitol	Amantadine	Amantadine (hydrochlorid e)	Clindamycin (hydrochlorid e)	Empty
g	Empty	Fenbendazole	Lincomycin (hydrochlorid e)	Moroxydine (hydrochlorid e)	Mycophenolic acid	Novobiocin	Novobiocin (sodium)	Quinine (hydrochlorid e dihydrate)	Ribavirin	Roxithromycin	Spectinomycin (dihydrochlorid)	Empty
h	Empty	Sulfadoxine	Tobramycin	Tobramycin (sulfate)	Phenytoin	Phenytoin (sodium)	Methacycline (hydrochlorid e)	Ciclopirox	Ciclopirox (olamine)	Econazole (nitrate)	Miconazole	Empty

Plate layout: CS-L002-20

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Miconazole (nitrate)	Lomefloxacin (hydrochlorid e)	Riboflavin	Cefprozil (monohydrate )	Azelastine (hydrochlorid e)	Oxacillin (sodium monohydrate)	Cloxacillin (sodium monohydrate)	Cloxacillin (sodium)	Amoxicillin (sodium)	Amoxicillin	Empty
b	Empty	Amoxicillin (trihydrate)	Neomycin (sulfate)	Streptomycin (sulfate)	Tetracycline (hydrochlorid e)	Thiamphenicol	Lonidamine	Clorsulon	Hygromycin B	Niclosamide	Niclosamide (olamine)	Empty
c	Empty	Danofloxacin (mesylate)	Enrofloxacin	Enrofloxacin (monohydrochlorid)	Nadifloxacin	Sulfathiazole	Sulfathiazole (sodium)	Ornidazole	Amikacin (hydrate)	Amikacin (disulfate)	Trimethoprim	Empty
d	Empty	Trimethoprim (lactate)	Sulfamerazine	Naftifine (hydrochlorid e)	Tylosin (tartrate)	Tylosin	Tylosin (phosphate)	Ampicillin	Ampicillin (sodium)	Ampicillin (trihydrate)	Flumequine	Empty

e	Empty	Azlocillin (sodium salt)	Ethambutol (dihydrochloride)	Clinafloxacin	Pentamidine (isethionate)	Probenecid	Norethindrone	Nafcillin (sodium monohydrate)	Ronidazole	Deferiprone	Zinc Pyrithione	Empty
f	Empty	Propranolol (hydrochloride)	Propranolol	Sulfacetamide (Sodium)	Cyclosporin A	Ceftazidime	Ceftazidime (pentahydrate)	Mafenide (Acetate)	Cefepime (Dihydrochloride)	Triclabendazole	Hyaluronic acid	Empty
g	Empty	Dirithromycin	Rabeprazole	Rabeprazole (sodium)	Vancomycin	Dapsone	Indinavir (sulfate)	Ecabet (sodium)	Cefepime	Ceftibuten (dihydrate)	Flomoxef (sodium)	Empty
h	Empty	Ceftriaxone	Ceftriaxone (sodium hydrate)	Pentoxifylline	Pazufloxacin (mesylate)	Pamidronate (disodium pentahydrate)	Itopride (hydrochloride)	Sertaconazole (nitrate)	Eflornithine	Eflornithine (hydrochloride hydrate)	L-Eflornithine (monohydrochloride)	Empty

Plate layout: CS-L002-21

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Fumagillin	Bicyclol	Artemotil	Cefozopran (hydrochloride)	Moxidectin	Milbemycin oxime	Acetazolamide	Acetazolamide (sodium)	Terbutaline (sulfate)	Terbutaline	Empty
b	Empty	Lumefantrine	Proguanil	Pivmecillinam (hydrochloride)	Pyridaben	Nitenpyram	Pymetrozine	Spirodiclofen	Dinotefuran	λ-Cyhalothrin	Chlorfenapyr	Empty
c	Empty	Thiophanate-Methyl	Metalaxyl	Prochloraz	Dimethomorph	Propiconazole	Tricyclazole	Azoxystrobin	Difenoconazole	Tebuconazole	Paclitaxel	Empty
d	Empty	Validamycin A	Cefmenoxime (hydrochloride)	Suramin (sodium salt)	Econazole	Permethrin	Ethacridine (lactate monohydrate)	Triacetin	Ceftiofur (sodium)	Bromperidol	Bromperidol (hydrochloride)	Empty
e	Empty	Rufloxacin hydrochloride	Tilmicosin	Tilmicosin (phosphate)	10-Undecenoic acid	Orbifloxacin	Propoxur	Azaserine	Succinylsulfathiazole	Danthron	Oxacillin (sodium salt)	Empty
f	Empty	Ethylparaben	Benzyl benzoate	Amprolium (hydrochloride)	Benzethonium (chloride)	Pidotimod	Nitromide	Sulfamonomoxine	Sulfantran	Oxethazaine	Paromomycin (sulfate)	Empty
g	Empty	Erythromycin Ethylsuccinate	Mupirocin	Chloramine-T	Sulfabenzamide	Cloxiquine	Thioridazine (hydrochloride)	Phthalylsulfacetamide	Cinchophen	Methicillin (sodium salt)	Dicloxacillin (Sodium hydrate)	Empty
h	Empty	Diethyltoluamide	Hexylresorcinol	Deferoxamine (mesylate)	Thiostrepton	Nithiamide	Cyanoacetohydrazide	Hexetidine	Chlorindanol	Oxolinic acid	Dinitolmide	Empty

Plate layout: CS-L002-22

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	8-Hydroxyquinoline	Edoxudine	Lanatoside C	Ftaxilide	Pasiniazid	Gemifloxacin (mesylate)	Flumethasone	Procodazole	Methylprednisolone succinate	Clindamycin phosphate	Empty
b	Empty	Lasalocid	Lasalocid (sodium)	Phenothrin	Morantel (tartrate)	Fosfomycin (calcium)	Cefazolin (sodium)	Tilorone (dihydrochloride)	Cinoxacin	Fenchlorphos	Aklomide	Empty
c	Empty	Hycanthone	Amitraz	Diazinon	Cefoxitin (sodium)	Secnidazole	Triclosan	Temephos	L-Cycloserine	Auranofin	Betamipron	Empty
d	Empty	Cefamandole (sodium)	Sulbentine	Imazalil	Broxaldine	Chlormidazole (hydrochloride)	Chlorhexidine (dihydrochloride)	Dioxanide furoate	Furaltadone (hydrochloride)	Bacampicillin (hydrochloride)	Clofoctol	Empty
e	Empty	Climbazole	Cephadrine	Imidazolidinyl urea	Nitroxoline	Bitoscanate	Lactulose	(+)-Camphor	Bekanamycin	(2S,5R,6R)-Ticarillin (disodium)	Crotamiton	Empty
f	Empty	Dixanthogen	Cefadroxil	Cefadroxil (hydrate)	Tetramisole (hydrochloride)	Piperonyl butoxide	Urethane	Pipemic acid	Dehydroacetic acid	Broxyquinoline	Trimipramine (maleate)	Empty
g	Empty	Bronopol	Sulfaphenazole	Flufenamic acid	Sisomicin (sulfate)	Exalamide	Ribostamycin (sulfate)	Acetohydroxamic acid	Dihydrostreptomycin (sulfate)	Dimetridazole	Thonzonium (bromide)	Empty
h	Empty	Chlorhexidine	Cefalonium	Cefuroxime (sodium)	Cefuroxime	Cefmetazole (sodium)	Chlorobutanol	Sulfaguanidine	Docosate (Sodium)	Cephalothin (sodium)	Cephalothin	Empty

Plate layout: CS-L002-23

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Sulfaquinoxaline	Sulfaquinoxaline (sodium salt)	Dimercaprol	Piperacillin (sodium)	Cetylpyridinium (chloride monohydrate)	2-Phenylethanol	Lithium citrate (tetrahydrate)	Ceforanide	Cefonicid (sodium)	Quinidine hydrochloride monohydrate	Empty
b	Empty	Carbaryl	Foscarnet (sodium)	Amodiaquine (dihydrochloride dihydrate)	Amodiaquine	Amodiaquine (dihydrochloride)	Oxiconazole nitrate	Cefuroxime axetil	Chlortetracycline (hydrochloride)	Apramycin (sulfate)	Cyromazine	Empty
c	Empty	Furazolidone	Carbadox	Oxantel (pamoate)	Piroctone olamine	Fusidic acid	Fusidic acid (sodium salt)	Oxyphenbutazone	Lincomycin (hydrochloride)	Methylene blue (trihydrate)	Chlorquinaldo	Empty



d	Empty	Meclocycline (Sulfosalicylate Salt)	Carbenoxolone (disodium)	Hydroxychloroquine (sulfate)	Florfenicol	Cefixime	Sulfamethoxy pyridazine	Lactitol (monohydrate)	Saccharin (sodium hydrate)	Diiodohydroxyquinoline	Phthalylsulfathiazole	Empty
e	Empty	Salicylanilide	Chloroxylenol	Tazobactam	Butylparaben	7-Aminocephalosporanic	Nifuroxazide	Acetarsol	Isoconazole (nitrate)	Clindamycin palmitate (hydrochloride)	Clindamycin	Empty
f	Empty	Sulconazole (mononitrate)	Cetylpyridinium (chloride)	Mezlocillin (sodium)	Ethoxzolamide	Moxalactam (sodium salt)	Alovudine	Thiacetazone	Citric acid (triammonium)	Benznidazole	(Z)-Capsaicin	Empty
g	Empty	Captan	Ceftizoxime	Cetalkonium (chloride)	Chloramphenicol palmitate	Sodium citrate (dihydrate)	Debrisoquin (hemisulfate)	Ditiocarb sodium	Ethyl Vanillate	Methenamine (hippurate)	Nifursol	Empty
h	Empty	Phenoxyethanol	Puromycin (dihydrochloride)	Pyridoxal phosphate	Quinidine (15% dihydroquinide)	Santonin	Spermine	Sulfachloropyridazine	Sulfamoxole	Sulfisomidin	Terconazole	Empty

Plate layout: CS-L002-24

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Thiethylperazine	Thiethylperazine (dimaleate)	Tosufloxacin (tosylate hydrate)	Triclocarban	Cefoxitin	Adefovir	Spectinomycin (dihydrochloride)	Oritavancin (diphosphate)	Fluazinam	Fuberidazole	Empty
b	Empty	Hexythiazox	Isoprothiolane	Kasugamycin (hydrochloride)	Kasugamycin (hydrochloride hydrate)	Metalddehyde	Faltan	Pirimiphosmethyl	Fenitrothion	Bromfenac (sodium)	Bromfenac (sodium hydrate)	Empty
c	Empty	(±)-Catechin	Cefazolin	Cefetamet pivoxil (hydrochloride)	Piperazine (phosphate)	Piperazine (tetraphosphate)	Methylprednisolone succinate	Diaveridine	Streptomycin	Rifamycin (sodium)	Midecamycin	Empty
d	Empty	Micronomicin	Acetylspiramycin	Josamycin	Piperacillin	Norvancomycin (hydrochloride)	Diniconazole	Thiacloprid	Iprodione	Thifluzamide	Fenvalerate	Empty
e	Empty	Famoxadone	Flusilazole	Carbosulfan	Phosalone	Pyriproxyfen	Pyrimethanil	Triflumuron	Tolclofomethyl	Tiamulin (fumarate)	Carboxin	Empty
f	Empty	Clofentezine	Cymoxanil	Cephaloridine (hydrate)	Etidocaine (hydrochloride)	Lucanthone	Sulfogaicol	Osalmid	Valpromide	Lactose	Ethopabate	Empty
g	Empty	Multifungin	Myclobutanil	Robenidine hydrochloride	Octenidine (dihydrochloride)	Carubicin (hydrochloride)	Ethacridine (lactate)	Sofalcone	Piperazine adipate	Thiamine nitrate	Thiamine disulfide	Empty
h	Empty	Sodium copper chlorophyllin	Lactate (calcium)	Hinokitiol	Povidone iodine	Lysozyme from chicken egg white	4-Methylherniarin	Quinine	2'-Deoxyuridine	Carbocysteine	ANS (ammonium)	Empty

Plate layout: CS-L002-25

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Malachite Green Carbinol base	Omeprazole metabolite Omeprazole	N-Desmethylclozapine	ITMN 4077	N-Acetylneuraminic acid	4-Aminosalicylic acid	Sodium 4-aminosalicylate (dihydrate)	2'-Aminoacetophenone	Picolinic acid	(-)-Epicatechin gallate	Empty
b	Empty	Honokiol	Oridonin	Curcumin	Demethoxycurcumin	Geniposide	Glycitin	Astragalol	Bergenol	Daidzin	Verbascoside	Empty
c	Empty	Arctigenin	Ginsenoside Rb1	Ginsenoside Rb2	Sanguinarin (chloride)	Psoralen	Osthole	Chlorogenic acid	Isochlorogenic acid A	3,4-Dicaffeoylquinic acid	4,5-Dicaffeoylquinic acid	Empty
d	Empty	Punicalagin	Solasodine	Ginkgolic Acid (C13:0)	(±)-Praeruptorin A	Betulinaldehyde	N6-Methyladenosine	Ancitabine (hydrochloride)	Guanosine	Isoliquiritigenin	Sophocarpine	Empty
e	Empty	Sophocarpine (monohydrate)	Palmitine (chloride)	Palmitine (hydroxide)	Dihydromyricetin	Hordenine	Xanthone	Sclareolide	Monensin sodium salt	Oleanolic Acid	Oxymatrine	Empty
f	Empty	Magnolol	Gramine	Cinchonidine	Dihydroartemisinin	Glycyrrhizic acid	Dipotassium glycyrrhizinate	Andrographolide	Artesunate	Asiatic acid	Azomycin	Empty
g	Empty	Baicalein	Baicalin	Ardisiacrispin A	Patchouli alcohol	Benzoic acid	Avicularin	Epigotritin	Bavachalcone	Psoralidin	Corylin	Empty
h	Empty	Attractylodin	Theaflavin	Saikosaponin A	Saikosaponin B2	Saikosaponin D	Hederacoside C	Cordycepin	Daphnetin	Isoimperatorin	Lycorine	Empty

Plate layout: CS-L002-26

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Lycorine (hydrochloride)	Paeoniflorin	Protocatechu aldehyde	Hederasaponin B	Octyl acetate	Allicin	Cholesterol	DL-Methionine	alpha-Mangostin	(+)-Magnoflorine (iodide)	Empty
b	Empty	3-Butylidenephthalide	Eugenol	Syringic acid	(E)-Ethyl p-methoxycinnamate	Methyl Paraben	Tuberostemone	Anthraquinone	Dihydrotanshinone I	(+)-Columbianetin	(+)-Columbianetin acetate	Empty

c	Empty	Falcarindiol	Sennoside A	Maackiain	Columbin	Glabridin	Harpagide	Artemether	Picroside II	Trigonelline (chloride)	Engeletin	Empty
d	Empty	Germacrone	2-Hydroxy-4-methoxybenzaldehyde	8-Gingerol	Hyperoside	Hypericin	L-Chicoric Acid	Corilagin	D-(-)-Quinic acid	Olaquinox	L-Lysine	Empty
e	Empty	L-Lysine hydrochloride	Phillyrin	α-Lipoic Acid	Gentiopicroside	Nitidine (chloride)	Erianin	Ethyl gallate	Pentagalloylgucose	(+)-Magnoflorine (chloride)	Xylitol	Empty
f	Empty	Cynaroside	(-)-α-Pinene	β-Pinene	Escin IA	Isoescin IA	Oroxilin A	Doxycycline	Doxycycline (hydrate)	Hydroxytyrosol	Purpurin	Empty
g	Empty	Pinocembrin	Solanesol	Deapioplatyodin D	Dehydrodiisoeugenol	Demethyleneberberine (chloride)	trans-Cinnamic acid	Mulberroside C	Morusin	Sorbic acid	Potassium sorbate	Empty
h	Empty	Maslinic acid	Eriodictyol	Dendrobine	Punicalin	α-L-Rhamnose monohydrate	L-Serine	Corypalmine	D-Pinitol	Usnic acid	(+)-Usnic acid	Empty

**Plate layout: CS-L002-27**

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Amentoflavone	Aucubin	Pseudolaric Acid C	Pseudolaric Acid A	Dehydrocorydalone	Dehydrocorydalone (chloride)	Dehydrocorydalone (hydroxyl)	Dehydroandrographolide	Dehydroandrographolide succinate	Kalidihydrographolide	Empty
b	Empty	Thiamine (hydrochloride)	α-Vitamin E	Schisandrin C	Schisandrin A	Daphnoretin	Agrimol B	Vanillic acid	Coumarin	Carvacrol	Berberine (sulfate)	Empty
c	Empty	Genkwanin	Coptisine (chloride)	Jatrorrhizine (chloride)	Jatrorrhizine (hydroxide)	Scutellarin	Scutellarein	(-)-Bornyl acetate	8-?O-Acetylharpagide	Isoferulic acid	Angelicin	Empty
d	Empty	Isoliquiritin	Isopimpinellin	Isomangiferin	Isolantalactone	Protopine	Tenuigenin	Narirutin	Camphor	Catalpol	Shikonin	Empty
e	Empty	Cephalotaxine	Bruceine A	Bevirimat	Dictamine	Alisol F	Gomisin G	Harringtonine	Catechin	Cimiracemide D	Dihydrosanguinarine	Empty
f	Empty	Dihydrochelythrine	Ajugol	Corydaline	(±)-Stylopine	(±)-Stylopine (hydrochloride)	Galegine (hydrochloride)	Coixol	beta-Mangostin	Epoxyzadira-dione	11-Ketofusidic acid	Empty
g	Empty	Hinokinin	Bleomycin A5	Xanthoxylin	Xanthohumol	Kaempferol-3-O-(6"-galloyl)-β-	Vasicine (hydrochloride)	Tricin	Mycaminosyllionolide	Thymidine	Strictosamide	Empty
h	Empty	Squalene	Sinapaldehyde	Robinetin	Isoscopoletin	Clinodiside A	(R)-Fangchinoline	Fangchinoline	Nevadensin	p-Anisic acid	20(R)-Ginsenoside Rh2	Empty

**Plate layout: CS-L002-28**

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	β-Caryophyllene	Pogostone	Citric acid	Citric acid (monohydrate)	Ferric citrate	Oxyresveratrol	Phytolaccagin	Methyl p-coumarate	Afzelin	Vincetoxicoside B	Empty
b	Empty	Betulonic acid	Hypocrellin B	Isoschaftoside	Picifeltarraenin IA	(-)-Fucose	Methyl palmitate	Oleanonic acid	Prunin	Pyrogallol	Halofuginone	Empty
c	Empty	Halofuginone (hydrobromide)	1-Methyl-2-pentyl-4(1H)-quinolinone	2,5-Dihydroxybenzaldehyde	2,6-Dimethoxy-1,4-	3,4-Dimethoxycinnamic acid	(E)-3,4-Dimethoxycinnamic acid	4'-O-Methylbavachalcone	Pristimerin	Isosinensetin	Nerolidol	Empty
d	Empty	cis-Nerolidol	Isoeugenol	5,7-Dihydroxychromone	(+)-(3R,8S)-Falcarindiol	Artemisic acid	Theaflavin 3,3'-digallate	Chebularic acid	Pyromeconic acid	Bellidifolin	Isoborneol	Empty
e	Empty	Ganoderic acid B	Methyl gallate	Octyl gallate	Verbenalin	Arteannuin B	Anacardic Acid	Propylparaben	Mosloflavone	(20R)-Protopanaxadiol	2,2':5',2"-Terthiophene	Empty
f	Empty	Evocarpine	Lindenenol	Cedrol	4',7'-Dimethoxyisoflavone	Apiopaeonoside	Polygalasaponin V	Deoxyshikonin	Swertianolin	4-Hydroxylonchocarpin	Isobavachromene	Empty
g	Empty	Picifeltarraenin IB	Guaijaverin	Rhapontigenin	N-p-trans-Coumaroyltyramine	Neuroscogenin	Cephaeline (dihydrochloride)	Gymnemagenin	Pleuromutilin	Piperlongumine	D-(+)-Melezitose	Empty
h	Empty	D-(+)-Melezitose (hydrate)	Procyanidin A2	Eudesmin	Chelidonine	Allosecurinine	Menthone	Febrifugine	Febrifugine (dihydrochloride)	Pinosylvin	Kukoamine A	Empty

**Plate layout: CS-L002-29**

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Stearyl glycyrrhettinate	Sinigrin (hydrate)	Paederoside	Tribuloside	Kakuol	(±)-Licarin A	Ganoderic acid G	Maltotetraose	Medicagenic acid	(E)-Methyl 4-coumarate	Empty

b	Empty	Lawsone	Deoxypodophyllotoxin	1-Monomyrustin	Diphyllin	Norwogonin	Rosamultin	(±)-Decursinol	Hypocrellin A	1-Kestose	Isoxanthohumol	Empty
c	Empty	Lupenone	Isoforsythiaside	Loureirin C	α-Lapachone	Sakuranetin	Bruceine D	Rupestonic acid	Artemitin	Grosvenorine	Tetrahydroepiberberine	Empty
d	Empty	Hexahydrofarnesyl acetone	Phytol	Ombuoside	Ombuin	Neotuberostemonine	Mulberrofuran G	Methyl 3-hydroxy-4,5-dimethoxyben	Matairesinol	Luteone	Loganetin	Empty
e	Empty	Lecanoric acid	Isodiospyrin	Canthin-6-one	Cauloside A	Coniferin	Curzerenone	Deoxylapachol	Fraxidin	Glaucine	Guaiol	Empty
f	Empty	Isoastilbin	Pseudolaric acid A-O-β-D-glucopyranoside	Trilobatin	5,7-Dihydroxy-4-methylcoumar	Dihydroartemisinic acid	Neodiosmin	17-Hydroxyisolatryhol	Lonicerin	3-O-Methylgalangin	Licoflavone B	Empty
g	Empty	Isopropyl ferulate	Yadanzoliolide A	Dehydrocorydaline (nitrate)	Kuwanon G	Coniferyl alcohol	Laetanine	Scutellarein tetramethyl ether	Thermopsine	Nepodin	Oxyberberine	Empty
h	Empty	DL-Goitrin	Estragole	Meloside A	2-Hydroxy-1-methoxyanthraquinone	Nonacosane	α-Terpineol	Methyl caffeate	Chryso-splenol D	8-Deoxygartanin	Aloin (mixture of A&B)	Empty

Plate layout: CS-L002-30

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Bacopasaponin C	Clemastanin B	Gartanin	(-)?-Maackiain	EGCG Octaacetate	Aflatoxin B1	Skullcapflavone II	Chlorothalonil	Pyraclostrobin	trans-Nerolidol	Empty
b	Empty	Luteolin-7-rutinoside	6-Hydroxycoumarin	Verbenone	Cefquinome (sulfate)	Vidarabine (monohydrate)	Virginiamycin M1	Calcimycin	Valinomycin	Apicidin	Beauvericin	Empty
c	Empty	β-Zearalenol	Citrinin	Penicillic acid	Patulin	Myricin	Netropsin (dihydrochloride)	Elemicin	Thymol	Retusin	4-Hydroxycoumarin	Empty
d	Empty	Dehydroabietic acid	Abietic acid	Rabdosin	Continentalic acid	Zingibroside R1	Usaramine	Sennidin B	Pseudolaric acid B β-D-glucoside	Pseudolaric Acid B	Hederacolchiside A1	Empty
e	Empty	Geraniol	Lapachol	α-Spinasterol	Cepharanthine	Ascr#2	Ascr#5	(-)-Corynoxidine	Zerumbone	7-Prenyloxycoumarin	11-Oxomogroside IIa	Empty
f	Empty	8-Epidiosbulbin E acetate	Nonanoic acid	Lactobionic acid	Helicin	Nerol	Difloxacin (hydrochloride)	Revaprazan (hydrochloride)	Mupirocin (calcium hydrate)	Maduramicin (ammonium)	Silymarin	Empty
g	Empty	Lanolin	D-?Arabinose	Benzoyleneurea	Ceftazidime (sodium)	Sulbenicillin (disodium)	Cefpodoxime Proxetil	Ceftiofur	Lactitol	Fenchyl alcohol	Sultamicillin (tosylate)	Empty
h	Empty	Leucomycin	Chloramphenicol succinate (sodium)	Sultamicillin	Lawsone methyl ether	Penicillin G procaine (hydrate)	Erythromycin estolate	Sulfacetamide	Cinnamyl acetate	Cinnamylideneacetic acid	Diphenylamine (hydrochloride)	Empty

Plate layout: CS-L002-31

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Penicillin G benzathine (tetrahydrate)	Spiramycin I	Citronellyl acetate	Daphnin	Obtusifolol	(E)-β-Farnesene	N-Hydroxypipercolic acid	DIMBOA	Neamine	Psoralenoside	Empty
b	Empty	Isopsoralenoside	12-O-Methylcarnosic acid	4'-Hydroxy-2,4-dimethoxycha	Acetylalkannin	5-Desmethylsinensetin	Tectol	Terminolic acid	Isookanin	Picropodophyllone	3,29-O-Dibenzoyloxykaroundiol	Empty
c	Empty	D-Dimannuronic acid (sodium)	L-Triguluronic acid	L-Tetraguluronic acid	L-Diguluronic acid (disodium)	Matairesinolside	Octanal	Colutehydroquinone	Protoneogracillin	Pangelin	Rapanone	Empty
d	Empty	Reuterin	Malabaricone B	Sodium Camptothecin	Enfufungin	Cichoriin	13,21-Dihydroeurycomanone	Granilin	Tellimagrandin II	Andrograpandin	Betulin diacetate	Empty
e	Empty	Piperitone	3-Methoxyflavone	2'-Fucosylactose	Aviptadil (acetate)	Aprotinin	Pepstatin	Pepstatin (Trifluoroacetate)	Enfuvirtide	Z-FA-FMK	Gramicidin	Empty
f	Empty	HA Peptide (TFA)	PBP10 (TFA)	Cyclosporin H	Pam3CSK4 (TFA)	LL-37, human	H-Lys-Trp-Lys-OH	HIV-1 TAT (48-60)	LCMV gp33-41	LCMV gp33-41 (TFA)	Nucleoprotein (396-404)	Empty
g	Empty	Nucleoprotein (396-404) (TFA)	TAT (48-57)	Ceratotoxin A	Temporin A	SPR741 (acetate)	Murepavadin (TFA)	Hepatitis B Virus Core (128-140)	Lactoferrin (17-41) (acetate)	HSV-gB2 (498-505)	Cyclo(L-Phe-L-Pro)	Empty
h	Empty	Cyclo(L-Pro-L-Val)	Cyclo(L-Leu-L-Pro)	Aureobasidin A	Dusquetide	Cyclo(L-Trp-L-Trp)	Tat-beclin 1	IDR-1	LLO (91-99)	Apelin-12	MDP1 (acetate)	Empty

Plate layout: CS-L002-32

1	2	3	4	5	6	7	8	9	10	11	12
---	---	---	---	---	---	---	---	---	----	----	----

a	Empty	Bulevirtide	Nisin Z	HCV Peptide (35-44)	Stilbamidine (dihydrochloride)	MGB-BP-3	Diffucortolone valerate	Tosylchloramide sodium trihydrate	BAY-Y 3118	Picloxydine	Tromantadine (hydrochloride)	Empty
b	Empty	PPA-904	Sulfabrom	Deprodone propionate	Furaprofen	Saperconazole	BF738735	Ethyl Orsellinate	1,3-Dithiane	3,5-Dimethoxybenzoic acid	4'-Hydroxy-3'-methylacetophenone	Empty
c	Empty	Glurate	5-Nitrobarbituric acid	2-Hydroxyacetophenone	Fluoroquinolonic acid	Rabepazole Sulfide	DHFR-IN-3	4-Methoxyphenethyl alcohol	Undecanoic acid	1-Heptadecanol	Gallic aldehyde	Empty
d	Empty	Phenazine (methylsulfate)	Abametapir	5-Nitro-1,10-phenanthroline	1-Hydroxy-2-butanone	Skatole	Methyl 2-amino-5-bromobenzoate	2-Aminoquinoline	Acridone	Phloracetophenone	Norharmane	Empty
e	Empty	2'-Deoxy-5'-O-DMT-2'-fluorouridine	N-(3-Oxoheptanoyl)-L-	Antimalarial agent 1	Hydroxycitric acid (tripotassium)	Tazobactam (sodium)	Bz-RS-iSer(3-Ph)-OMe	Diazolidinyl urea	1-Deoxymannojirimycin	3,4,5-Trimethoxybenzaldehyde	Methylisothiazolinone (hydrochloride)	Empty
f	Empty	Benzylurea	4'-Hydroxy-2'-methylacetophenone	Methylisothiazolinone	trans-4-Methylcyclohexanamine	2-Chloroacetamide	Isoxazole	Quinine (sulfate hydrate)	Fmoc-N-Me-Phe-OH	Dodecylmethylpentasiloxane	Cyclofenil	Empty
g	Empty	2'-Deoxy-2'-fluoroguanosine	Taurolidine	Tulobuterol (hydrochloride)	2'-O-Methylcytidine	Muramic acid	2'-Deoxy-2'-fluorocytidine	2,6-Dichlorodiphenylamine	3'-Deoxyuridine	2,2'-Anhydruvidine	2-Hydroxycinnamic acid	Empty
h	Empty	Decyl aldehyde	Benzothiazole	DL-Pyrogutamic acid	3-Nitropropanoic acid	Trimethyloctadecylammonium (bromide)	Foscarnet (trisodium hexahydrate)	N4-Acetylsulfamethoxazole	3-Deazaadenosine	2'-Deoxy-2'-fluorouridine	2',3'-Dideoxyadenosine	Empty

Plate layout: CS-L002-33

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	4-Isopropylbenzoic acid	4-Methoxycinnamaldehyde	Chlorhexidine (diacetate)	Thianthrene	Eugenol acetate	Isobutylparaben	Desaminotyrosine	4-Methylcinnamic acid	Mandelic acid	4-Methyl-5-thiazolethanol	Empty
b	Empty	T-1105	2-Benzoxazolinone	Phthalide	2-Heptanol	2-Amino-2'-deoxyadenosine	L-Methioninamide	Naphthalene-2,3-Dicarboxaldehyde	Fosfomycin (sodium)	Isopropyl 4-hydroxybenzoate	Acivicin (hydrochloride)	Empty
c	Empty	5-Fluoroarotic acid	4-Chlorosalicylic acid	2-Undecanone	DL-3-Phenylactic acid	Methyl cinnamate	1-Naphthalene methanol	Tiliquinol	D-Cysteine	N-Phenylacrylamide	Bifendate	Empty
d	Empty	4(3H)-Quinazolinone	Albendazole sulfone	Flutriafol	Pyribencarb	7-Chloro-4-(piperazin-1-yl)quinoline	α-Terpinene	Triflumizole	Quipazine	Hydroxychloroquine	4-Chloroguaiacol	Empty
e	Empty	2,4-Dichlorobenzyl alcohol	Methyl α-D-mannopyranoside	Nifurtimox	Kanamycins (sulfate)	Chromomycin A3	N2-Acetylguanidine	3,5-Di-tert-butylphenol	Fmoc-Glu-OMe	Aegeline	N-Decyl-N,N-dimethyldecylamine-1-aminium	Empty
f	Empty	Oxychloroquine	N-Hexanoyl-DL-homoserine	Nitroxynil	Carvacrol methyl ether	Kojic acid	Antitrypanosomal agent 1	Etozazole	2-Aminoimidazole	5,7-Dihydroxycoumarin	Antibacterial agent 18	Empty
g	Empty	(S)-Tenofovir	3-Pentanol	Acriflavine (hydrochloride)	2-Ethyl-6-methylphenol	Heneicosane	Chlorobutanol (hemihydrate)	3'-O-Methylguanosine	4-Hydroxyhygric acid	5'-O-DMT-2'-TBDMS-Uridine	Ferron	Empty
h	Empty	Arprinocid	Transtorine	Etzadroxil	APOBEC3G-IN-1	Antibacterial agent 117	Isatin-β-thiosemicarbazone	Salinazid	HIV-1 inhibitor-48	BVDU 5'-Triphosphate (ammonium)	L-Norleucine	Empty

Plate layout: CS-L002-34

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	Thiosemicarbazide	5-Hydroxymethylfurfural	Phenothiazine	Oxindole	4-Hydroxyacetophenone	Cinchonine	Farnesol	Benzocaine	4-Hydroxybenzoic acid	Saccharin	Empty
b	Empty	Acetic acid (potassium 99%)	Acetic acid (magnesium tetrahydrate)	Lauric acid	Maleic Acid	Semicarbazide (hydrochloride)	L-Lactic acid	HODHBt	DL-Serine	Potassium Gluconate	trans-Chalcone	Empty
c	Empty	Pipecolic acid	4-Chlorocinnamic acid	(1R)-α-Pinene	Ethyl acetoacetate	Tridecanoic acid	2-Carboxybenzaldehyde	2-Phenylacetophenone	o-Vanillin	Empty	Empty	Empty
d	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty
e	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty
f	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty
g	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty
h	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty