


=====

Acq. Operator	: Su Xiao Ying(LCMS-02)	Seq. Line	: 41
Acq. Instrument	: HY-LCMS-02	Location	: P1-C-06
Injection Date	: 4/27/2016 1:55:47 PM	Inj	: 1
		Inj Volume	: 3.000 µl
Acq. Method	: D:\AGLIENT 1260\DATA\20160427\20160427 2016-04-27 09-26-16\100-1000MS+3MIN-1.5_(0.02%FA).M		
Last changed	: 4/27/2016 9:26:16 AM by Su Xiao Ying(LCMS-02)		
Analysis Method	: D:\AGLIENT 1260\DATA\20160427\20160427 2016-04-27 09-26-16\100-1000MS+3MIN-1.5_(0.02%FA).M (Sequence Method)		
Last changed	: 4/27/2016 2:06:23 PM by Su Xiao Ying(LCMS-02) (modified after loading)		
Method Info	: HY-365_5H01RS,M,A-RP-108, 210nm,23min		
Catalog No	: CS-0864 Batch#20227 A-RP-134		

Additional Info : Peak(s) manually integrated



=====
Area Percent Report
=====

Sorted By : Signal
Multiplier : 1.0000
Dilution : 1.0000
Do not use Multiplier & Dilution Factor with ISTDs

Signal 1: DAD1 B, Sig=214,4 Ref=off


Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	1.739	MM	0.0573	8161.13623	2372.68750	98.5073
2	2.563	MM	0.0679	123.66550	30.33958	1.4927

Totals : 8284.80173 2403.02708

=====
*** End of Report ***


```
=====
Acq. Operator   : Su Xiao Ying(LCMS-02)          Seq. Line :   41
Acq. Instrument : HY-LCMS-02                     Location  : P1-C-06
Injection Date  : 4/27/2016 1:55:47 PM           Inj       :    1
                                                Inj Volume: 3.000 µl
Acq. Method     : D:\AGLIENT 1260\DATA\20160427\20160427 2016-04-27 09-26-16\100-1000MS+3MIN-
                  1.5_(0.02%FA).M
Last changed    : 4/27/2016 9:26:16 AM by Su Xiao Ying(LCMS-02)
Analysis Method : D:\AGLIENT 1260\DATA\20160427\20160427 2016-04-27 09-26-16\100-1000MS+3MIN-
                  1.5_(0.02%FA).M (Sequence Method)
Last changed    : 4/27/2016 2:08:33 PM by Su Xiao Ying(LCMS-02)
                  (modified after loading)
Method Info     : HY-365_5H01RS,M,A-RP-108, 210nm,23min
Catalog No      : CS-0864 Batch#20227
                  A-RP-134
=====
```

Additional Info : Peak(s) manually integrated



MS Signal: MSD1 TIC, MS File, ES-API, Pos, Scan, Frag: 50
Spectra averaged over upper half of peaks.
Noise Cutoff: 1000 counts.
Reportable Ion Abundance: > 10%.

Retention Time (MS)	MS Area	Mol. Weight or Ion
1.748	5612561	308.20 I
		307.20 I



*** End of Report ***