

PROTAC Library (96-well)

Product Details:

Catalog Number:	CS-L151
Formulation:	A collection of 61 PROTACs supplied as pre-dissolved Solutions or Solid.
Container:	96- or 384-well Plate with Peelable Foil Seal; 96-well Format Sample Storage Tube With Screw Cap and Optional 2D Barcode
Storage:	In solvent: -80°C 2 years, -20°C 1 year; Powder: -20°C 3 years, 4°C 2 years
Shipping:	Blue ice
Packaging:	Inert gas

Plate layout: CS-L151-1

	1	2	3	4	5	6	7	8	9	10	11	12
a	Empty	ARV-771	dBET1	PROTAC CDK9 Degradere-1	PROTAC BET Degradere-1	MS4078	MZP-54	MZP-55	dBET6	PROTAC BET Degradere-10	PROTAC BET degradere-2	Empty
b	Empty	MD-224	SJF α	PROTAC HSP90 degradere BP3	dBRD9	TL13-12	CCT367766 (formic)	TL13-112	GMB-475	ERD-308	SIAIS178	Empty
c	Empty	PROTAC CRBN Degradere-1	XY028-133	KB02-SLF	ARCC-4	MS432	TL12-186	TD-165	FKBP12 PROTAC RC32	GSK215	ZXH-4-130 (TFA)	Empty
d	Empty	SJF620	dCBP-1	MD-222	BSJ-03-204	BSJ-04-132	CMP98	CRBN-6-5-5-VHL	PROTAC RAF degradere 1	MI-389	XY028-140	Empty
e	Empty	BSJ-4-116	SK-575	PROTAC ERR α Degradere-3	SIM1	JB170	PROTAC-O412	NJH-2-056	AU-15330	PROTAC AR-V7 degradere-1	dTAGV-1 (TFA)	Empty
f	Empty	JQAD1	(S,R,S)-AHPC-C2-amide-	dTAG-47	PROTAC EZH2 Degradere-1	MS9427	MS9427 (TFA)	MS177	MS8815	YX-2-107	MTX-23	Empty
g	Empty	ARV-825	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty
h	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty	Empty